ENGLISH

ZIMMER® DERMATOME II BLADE INSTRUCTION MANUAL

INDICATIONS FOR USE

The Zimmer Dermatome II Blade, when used with the Zimmer® Dermatome II, is intended to provide variable graft thickness and width capabilities.

DESCRIPTION

The Zimmer Dermatome II Blade is a stainless steel blade attached to a molded mount.

WARNINGS AND PRECAUTIONS

- For use only with the Zimmer Dermatome II.
- Warning! Blade is sharp. Handle with care.
- Reuse of a single use device that has come in contact with blood, bone, tissue or other body fluids may lead to patient or user injury. Possible risks associated with reuse of a single use device include, but are not limited to, mechanical failure and transmission of infectious agents.
- Sterile only if package is unopened and undamaged.

BLADE INSTALLATION

Lubrication of the blade is not necessary. Use a new, sterile blade for each procedure. Additional blades may be required for extensive procedures.

- Place the Zimmer Dermatome II handpiece throttle in the SAFE (O•) position. To place the handpiece in the SAFE position, slide the safety lock on the throttle toward the blade end of the instrument to the SAFE (O•) position. Only the SAFE symbol (O•) should be visible (see Figure 1).
- Using the Zimmer Dermatome II screwdriver (REF 00-8851-205-00), loosen the handpiece width plate screw. Do not remove the screw from the handpiece. If a width plate is already attached, remove by holding both sides and sliding/lifting forward and up (Figure 2). Do not lift from the front as this will cause contact with the blade and possible injury.
- If replacing a blade, remove the used blade before inserting a new one. Refer to BLADE REMOVAL section.
- Place entire blade/cover assembly (see Figure 3) onto the handpiece such that the blade key fits into the notch of the handpiece, and the handpiece drive pin mates with the hole in the blade (see Figure 4). Place your finger on the blade assembly (white portion). While holding the blade assembly, lift one end of the blade cover (blue portion) and remove blade cover from blade assembly (see Figure 5). **Retain blade cover. The blade cover can be used for blade removal from the device after use.**
- With the printing on the width plate facing out, place the width plate flat on the face of the dermatome head (see Figure 6). Carefully slide back toward the screw, until the width plate hole is located under the mounting screw. Do not grasp the front of the width plate as this will cause contact with the blade and possible injury. Using the screwdriver, gently tighten the mounting screw. Do not over-tighten screw (see Figure 7).
- NOTE: Choose the proper width plate to satisfy cutting requirements.

BLADE REMOVAL

- Using the screwdriver, loosen the width plate screw. Do not remove the screw.
- To remove the width plate, hold both sides of the width plate and slide forward and lift up. Do not lift from the front as this will cause contact with the blade and possible injury (Figure 2).
- Place the blade cover over one end of the blade assembly and snap it down onto other end of the blade assembly (preferred, see Figure 8), or remove the blade assembly carefully by lifting up from under the blade assembly using the finger access area (see Figure 9).
- Dispose of used blade in a sharps container.

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© 2011 Zimmer Surgical, Inc.
ATTENZIONE: la legge federale (U.S.A.) limita la vendita di questo dispositivo a medici o dietro prescrizione medica.
Zimmer Surgical, Inc.
Dover, Ohio 44622 U.S.A.
Printed in U.S.A
Zimmer U.K. LTD.